

PROCATION

0 O°

(V) ANO CAL

Endo it your way!

Find the VDW endo motor that matches your style

VDW.GOLD

ROCATION

0

VDW.CONNECT Connectable endo devices

That's how efficient endo works today!



VDW.CONNECT Drive® Smart endo motor

VDW.CONNECT® App Smart treatment monitoring VDW.CONNECT Locate® Smart apex locator

My ON

VDW.CONNECT Drive®

Simplicity-driven to make preparation smoother and safer

- Cordless and ergonomic handpiece for maximum comfort and mobility
- Easily set for a wide range of rotary and reciprocating file systems through the VDW.CONNECT[®] App*
 - * without app preset for preparation with RECIPROC[®] blue, RECIPROC[®] and R-PILOT[™]

360° rotating mini-head for excellent visibility

Apple, the Apple logo and iPad® are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc.. The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth® SIG, Inc..

VDW.CONNECT Locate®

Precision at your fingertips to improve your treatment

- Mini size, lightweight and connected via Bluetooth® to integrate easily in your workflow
- Clear LED and acoustic signals
- Easily control and visualize the file progression on the iPad[®] screen through the VDW.CONNECT[®] App
- Shaping target and apical reverse function (app required)



VDW.CONNECT® App

The intuitive way to better results

- Comprehensive support throughout the entire endo process with VDW.CONNECT Drive® and VDW.CONNECT Locate®
- Personalize your settings, customize file sequences and record treatments
- Get quick access to patient education, tips and tricks, DFU and latest updates
- Change from reciprocating to rotary




VDW.GOLD®RECIPROC®

Endo motor with integrated apex locator A market proven device you can rely on!

VDW.GOLD



- RECIPROC REVERSE Comfort Function**: two-stage indicator helps determine the right method of preparation
- Additional setting for root canals with difficult anatomy for rotary motion



ANA*: Reduction of torque limit during the preparation of canals with difficult anatomies

ASR*: Automatic reverse rotation when the set torque limit is reached

*in rotary mode only **in RECIPROC® mode only

High performance for rotary and reciprocating root canal preparation

- Powerful and safe preparation in torque controlled rotary and reciprocating motion
- Supports an extensive range of rotary and reciprocating file systems
- Dr's Choice: store up to 15 individual torque and speed settings
- Battery-operated motor; maintains constant operation even during charging

Powerful preparation and precise length determination in one device

Choose between three working modes:

- **1.** Monitor the file progression simultaneously during preparation.
- **2.** Prepare the root canal without length determination.
- **3.** Determinate the working length or monitor the file progression independently of preparation.

21	-		
V	0	V	

THE GOLD SUCCESS STORY

In 2007, VDW challenged the established endo standards by introducing VDW.GOLD[®] as the first endo motor with an integrated apex locator. In 2011, the scope was extended to the innovative preparation method RECIPROC[®] with the new VDW.GOLD[®]RECIPROC[®] edition. Today, VDW's GOLD motor is one of the best-selling endo motors on the market, always reflecting VDW's pure endo expertise and German heritage of engineering excellence. Continuous improvements make it one of the most reliable endo motors, coming with an exceptional 3-year guarantee. An outstanding success that VDW shares with thousands of VDW.GOLD[®] users all over the world.



VDW.SILVER®RECIPROC®

Stand-alone endo motor Essential performance for effective and safe root canal preparation

- Choose from a wide range of reciprocating and rotary file systems
- Dr's Choice: store up to 15 individual torque and speed settings
- ASR*: automatic reverse rotation when the set torque limit is reached
- RECIPROC REVERSE Comfort Function**: two-stage indicator identifies the right method of preparation
- Battery-operated motor; maintains constant operation even during charging

VDW.SILV

*in rotary mode only **in RECIPROC® mode only

RAYPEX®6

Stand-alone apex locator Length determination with user-friendly interface

- Clear acoustic and visual signals on a bright color touch screen
- Various safety functions to enhance user and patient confidence
- Demo mode for easy and comprehensible patient education
- Intuitive handling, handy pocket-sized design, foldable
- Long-lasting rechargeable battery

A powerful duo for professional endo treatments

VDW.SILVER®RECIPROC® & RAYPEX®G

RECIPROC[®] direct

The first reciprocating contra-angle It's not a big deal to start with RECIPROC[®]

- A quick and inexpensive way to the most time-saving preparation technique
- Simply click it on your electric or air-driven motor*
- Drives RECIPROC® instruments safely and efficiently
- Ideal accessory to your dental unit allowing you to switch easily from manual to mechanical preparation
 - * ISO-E coupling according to ISO 3964 is necessary

Be pragmatic: Simply make your dental unit an endo unit!

My choice? Plug'n'play! RECIPROC direct

Find the VDW endo motor that matches your style!

		VDW. CONNECT Drive [®] stand- alone	VDW. CONNECT Drive® + App + Locate®	VDW. GOLD® RECIPROC®	VDW. SILVER® RECIPROC®	RECIPROC [®] direct
with s rotray	brary specific pre-settings of reciprocating and NITi systems					
Reciprocating systems	RECIPROC [®] , RECIPROC [®] blue, R-PILOT™	\checkmark	\checkmark	\checkmark	\checkmark	\checkmark
Recipr syst	WaveOne®, WaveOne® Gold, WaveOne® Gold Glider	×	✓	✓	✓	×
Continous rotary NiTi systems	VDW.ROTATE [™] , VDW.ROTATE [™] Retreatment, Mtwo [®] , Mtwo [®] Retreatment, FlexMaster [®] , ProTaper [®] , ProTaper Next [®] , ProTaper [®] Universal, ProTaper [®] Retreatment, ProGlider [®]	×	~	~	~	×
ontinou NiTi sy:	K3™, Gates	×	×	\checkmark	✓	×
O	ProFile®, GT® Series X, Hero®, RaCe, FlexMaster® Retreatment, Lentulo	×	×	~	×	×
Indiv	idually programmable torque and speed settings	×	~	~	~	-
Indiv	idually programmable file sequences	×	~	×	×	-
Special safety featu (in rotary mode)	ASR: Automatic reverse rotation when the set torque limit is reached	×	~	✓	✓	×
Special se (in rot	ANA Programme: Reduction of torque limit enhances safety during the preparation of canals with difficult anatomies	×	~	~	×	×
Two-	PROC REVERSE Comfort Function stage indicator assists the right method of aration (in RECIPROC® mode)	~	~	~	~	×
Softv	vare update possible	\checkmark	\checkmark	✓	✓	-
Cons	tant operation even during charging	\checkmark	\checkmark	✓	✓	-
Meni Displ File s	u ay can be individually customized: ystems that are not in use can be deactivated.	-	~	~	×	-
• Cho sep • Aut	Iocator Dice between simultaneous length control and Derate length determination Dematic apex stop during preparation with Procating and rotary systems	×	~	~	×	×
6:1 c	ontra-angle	\checkmark	\checkmark	✓	✓	-



For an endo treatment in the original reciprocating motion, VDW recommends using RECIPROC[®], RECIPROC[®] *blue* and R-PILOT[™] files in combination with the VDW motors VDW.CONNECT Drive[®], VDW.GOLD[®]RECIPROC[®], VDW.SILVER[®]RECIPROC[®] and RECIPROC[®] direct.

WaveOne®, ProTaper®, ProTaper Next®, ProGlider® and GT® are registered trademarks owned by Dentsply Sirona. K3™ is a trademark owned by KAVO KERR GROUP. Hero® is a registered trademark owned by Micro Mega.

Developed to work together

Order information

■ VDW.CONNECT Drive®

VDW.CONNECT Drive [®] Incl. 6:1 contra-angle and handpiece holder	- V04 1183 000 000 (DFU printed in English) V04 1184 000 000 (DFU printed in TK, KR, CHI, BR, RUS, UKR, IND)
 VDW.CONNECT Drive[®] Set	- V04 1183 000 500
 VDW.CONNECT Drive[®] Set Medium	_ V04 1181 000 000
 VDW.CONNECT Drive[®] Set Medium blue	

- User card, brochure, directions for use

Accessories for VDW.CONNECT Drive[®]

VDW.CONNECT Drive [®] Handpiece without contra-angle	V04 1185 000 502
6:1 contra-angle	V04 1185 000 504
Handpiece sleeves, 100 pieces	V04 1185 000 503
Charger incl. adapters	V04 1185 000 501
Spray nozzle for manual lubrication	V04 1185 000 506
Contra-angle adapter for Assistina	V04 1185 000 511 (corresponds with W&H REF 06037100)

VDW.CONNECT Locate®

VDW.CONNECT Locate® Kit	t_
-------------------------	----

• 1 VDW.CONNECT Locate® Apex Locator

- Disposable barrier sleeves
- Clamp, measurement cable, external tester
- $2 \times \text{file clip}, 5 \times \text{lip clip}$
- Adapter for VDW.CONNECT Drive® Handpiece
- Rechargable battery and charger
- Directions for use, Getting Started Guide

M00 CLOC 000 000

Order information

■ VDW.CONNECT Drive[®] and Locate[®]

 VDW.CONNECT Drive® & Locate® Kit VDW.CONNECT Drive incl. 6:1 contra-angle and handpiece holder VDW.CONNECT Locate® Kit 	M00 CLOC CDR KIT
 VDW.CONNECT Drive® & Locate® RECIPROC® blue VDW.CONNECT Drive incl. 6:1 contra-angle and handpiece holder VDW.CONNECT Locate® Kit RECIPROC® blue System Kit 24: 24 RECIPROC® blue instruments (25 mm), sterile, 12 × R25, 6 × R40, 6 × R50 180 RECIPROC® blue Paper Points, assorted, sterile, 60 × R25, 60 × R40, 60 × R50 60 RECIPROC® blue Gutta-Percha cones, assorted, 36 × R25, 12 × R40, 12 × R50 1 Endo training tooth, trepanated 	M00 CLOC RBL KIT

- User card, brochure, directions for use

Accessories for iPad[®]

VDW.CONNECT iPad® Stand	M00 CO00 STA 000
Barrier sleeves for iPad mini [®] 4, 100 pieces	M00 CO00 079 BAR
Barrier sleeves for iPad [®] 9,7", 100 pieces	M00 CO00 097 BAR
Barrier sleeves for iPad [®] 10,5", 100 pieces	_ M00 CO00 105 BAR
Barrier sleeves for iPad [®] 12,9", 100 pieces	M00 CO00 129 BAR

■ VDW.GOLD®RECIPROC®

VDW.GOLD®RECIPROC® Incl. 6:1 contra-angle	V04 1173 000 000
 VDW.GOLD®RECIPROC® + RECIPROC® System Kit 12 RECIPROC® instruments (25 mm), sterile, 6 × R25, 3 × R40, 3 × R50 144 RECIPROC® Paper Points, assorted, sterile, 96 × R25, 24 × R40, 24 × R50 60 RECIPROC® Gutta-Percha cones, assorted, 40 × R25, 10 × R40, 10 × R50 User card, brochure, directions for use 	V04 1173 025 611
 VDW.GOLD®RECIPROC® <i>blue</i> 24 VDW.GOLD®RECIPROC® RECIPROC® <i>blue</i> System Kit 24: 24 RECIPROC® <i>blue</i> instruments (25 mm), sterile, 12 × R25, 6 × R40, 6 × R50 180 RECIPROC® <i>blue</i> Paper Points, assorted, sterile, 60 × R25, 60 × R40, 60 × R5 60 RECIPROC® <i>blue</i> Gutta-Percha cones, assorted, 36 × R25, 12 × R40, 12 × R50 1 Endo training tooth, trepanated User card, brochure, directions for use 	0
VDW.SILVER®RECIPROC®	
VDW.SILVER®RECIPROC® Incl. 6:1 contra-angle	V04 1163 000 000
VDW.SILVER®RECIPROC® + RECIPROC® System Kit 12 RECIPROC® instruments (25 mm), sterile, 6 × R25, 3 × R40, 3 × R50	V04 1163 025 611

- 144 RECIPROC® Paper Points, assorted, sterile, 96 × R25, 24 × R40, 24 × R50
- 60 RECIPROC® Gutta-Percha cones, assorted, 40 × R25, 10 × R40, 10 × R50
- 1 Endo training tooth, trepanated
- User card, brochure, directions for use

Accessories for VDW.GOLD®RECIPROC® and VDW.SILVER®RECIPROC®

6:1 Contra-angle	_ V04 1079 000 000
Contra-angle sleeve, 5 pieces	_ V04 1107 000 512
Lip clip, 5 pieces	_ V04 1107 000 514
File clip, 2 pieces	V04 1107 000 515
File clip cable	_ V04 1177 000 516
Lip clip cable with ferrite ring	_ V04 1107 000 519
Charger	_ V04 1107 000 517
Foot switch	_ V04 1107 000 510

■ RAYPEX®6

RAYPEX®6	V04 1113 000 000
Incl. accessories	

■ RECIPROC[®] direct

RECIPROC [®] direct <i>blue</i>	21 mm V04 1094 021 624*
 1 RECIPROC[®] direct contra-angle 24 RECIPROC[®] blue instruments, sterile, 12 × R25, 6 × R40, 6 × R50 1 spray nozzle for manual lubrication Directions for use, user card *only for EOC and Japanese markets 	25 mm V04 1095 025 000 V04 1094 025 624*
RECIPROC [®] direct Starter Kit** 1 RECIPROC [®] direct contra-angle	21 mm V04 1093 021 000
 44 RECIPROC® instruments, sterile, 32 × R25, 8 × R40, 4 × R50 1 spray nozzle for manual lubrication Directions for use, user card **not available in Europe 	25 mm V04 1093 025 000
 RECIPROC[®] direct Contra-Angle	V04 1091 000 000

• Directions for use, user card



For further information please contact your preferred distributor or visit www.vdw-dental.com.

WDW GmbH

Bayerwaldstr. 15 81737 Munich Germany Tel +49 (0) 89 62734-0 Fax +49 (0) 89 62734-304 info@vdw-dental.com

